

Product datasheet for **MR224415L4V**

IL5ra (NM_008370) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	IL5ra (NM_008370) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	IL5ra
Synonyms:	CD125; CDw125; IL5r
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_008370
ORF Size:	1245 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224415).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008370.2 , NP_032396.1
RefSeq Size:	3557 bp
RefSeq ORF:	1248 bp
Locus ID:	16192
UniProt ID:	P21183
Cytogenetics:	6 49.19 cM
Gene Summary:	This is the receptor for interleukin-5. The alpha chain binds to IL5.[UniProtKB/Swiss-Prot Function]



[View online »](#)