

## Product datasheet for **MR224428L4V**

### Abcc3 (NM\_029600) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Abcc3 (NM_029600) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Abcc3
Synonyms:	1700019L09Rik; ABC31; MLP2; MOAT-D; MRP3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_029600
ORF Size:	4566 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224428).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_029600.3</a>
RefSeq Size:	5002 bp
RefSeq ORF:	4569 bp
Locus ID:	76408
UniProt ID:	<a href="#">B2RX12</a>
Cytogenetics:	11 D
Gene Summary:	May act as an inducible transporter in the biliary and intestinal excretion of organic anions. Acts as an alternative route for the export of bile acids and glucuronides from cholestatic hepatocytes.[UniProtKB/Swiss-Prot Function]



[View online »](#)