

Product datasheet for **MR224453L2V**

Lgals3bp (NM_011150) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lgals3bp (NM_011150) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lgals3bp
Synonyms:	90K; CyCAP; MAC-2BP; Ppicap; Tango10b
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_011150
ORF Size:	1731 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224453).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011150.2
RefSeq Size:	2168 bp
RefSeq ORF:	1734 bp
Locus ID:	19039
UniProt ID:	Q07797
Cytogenetics:	11 E2
Gene Summary:	Promotes integrin-mediated cell adhesion. May stimulate host defense against viruses and tumor cells (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)