

Product datasheet for **MR224459L4V**

Phf20 (NM_172674) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Phf20 (NM_172674) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Phf20
Synonyms:	6820402O20Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_172674
ORF Size:	3030 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224459).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_172674.2 , NP_766262.2
RefSeq Size:	5752 bp
RefSeq ORF:	3033 bp
Locus ID:	228829
UniProt ID:	Q8BLG0
Cytogenetics:	2 H1



[View online »](#)

Gene Summary:

Contributes to methyllysine-dependent p53/TP53 stabilization and up-regulation after DNA damage (By similarity). Methyllysine-binding protein, component of the MOF histone acetyltransferase protein complex. Not required for maintaining the global histone H4 'Lys-16' acetylation (H4K16ac) levels or locus specific histone acetylation, but instead works downstream in transcriptional regulation of MOF target genes. As part of the NSL complex it may be involved in acetylation of nucleosomal histone H4 on several lysine residues. [UniProtKB/Swiss-Prot Function]