

Product datasheet for MR224470L4

Jazf1 (NM_173406) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Jazf1 (NM_173406) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Jazf1
Synonyms:	AI591476; C820002C15; Jaz1; Tip27
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224470).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_173406
ORF Size:	729 bp



[View online »](#)

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	<p>This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.</p>
Components:	<p>The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).</p>
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<p>NM_173406.3, NP_775582.2</p>
RefSeq Size:	<p>2925 bp</p>
RefSeq ORF:	<p>732 bp</p>
Locus ID:	<p>231986</p>
UniProt ID:	<p>Q80ZQ5</p>
Cytogenetics:	<p>6 25.74 cM</p>
Gene Summary:	<p>Acts as a transcriptional corepressor of orphan nuclear receptor NR2C2 (By similarity). Inhibits expression of the gluconeogenesis enzyme PCK2 through inhibition of NR2C2 activity (PubMed:24380856). Also involved in transcriptional activation of NAMPT by promoting expression of PPARA and PPARC (PubMed:24930994). Plays a role in lipid metabolism by suppressing lipogenesis, increasing lipolysis and decreasing lipid accumulation in adipose tissue (PubMed:24380856, PubMed:25614086). Plays a role in glucose homeostasis by improving glucose metabolism and insulin sensitivity (PubMed:25614086, PubMed:24380856). [UniProtKB/Swiss-Prot Function]</p>

Product images:



Circular map for MR224470L4