

## Product datasheet for MR224475

### Ifnl2 (NM\_001024673) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ifnl2 (NM\_001024673) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Ifnl2  
**Synonyms:** EG330496; IL-28A; IL28a  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR224475 representing NM\_001024673  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCCGGATCGCC

ATGCTCCTCTGCTGTTGCCTCTGCTGCTGGCCGAGTGTGACAAGAACCCAAGCTGACCCGTGCCCA  
 GGGCCACCAGGCTCCCAGTGAAGCAAAGGATTGCCACATTGCTCAGTTCAAGTCTGTGCCCAAAAGA  
 GCTGCAGGCTTCAAAAAGGCCAAGGATGCCATCGAGAAGAGGCTGCTTGAGAAGGACCTGAGGTGCAGT  
 TCCCACCTTTCAGGCTGGGACCTGAAGCAGCTGCAGGTCCAAGAGCGCCCAAGGCCTTGACAG  
 CTGAGGTGGCCCTGACCCTGAAGGTCTGGGAGAACATGACTGACTCAGCCCTGGCCACCATCCTGGCCA  
 GCCTTTCATACACTGAGCCACATTCACTCCCAGCTGCAGACCTGTACACAGTTTCAAGGCCACAGCAGAG  
 CCCAGGTCCCGAGCCGCGCCTCTCCCGCTGGCTGCACAGGCTCCAGGAGGCCAAGGAGACCC  
 CTGGCTGCCTGGAGGCTCTGTACCTCCAACCTGTTTCGCTGCTACCCGGGACCTCAAGTGTGTGGC  
 CAATGGAGACCAAGTGTGTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR224475 representing NM\_001024673  
 Red=Cloning site Green=Tags(s)

MLLLLLPLLLA AVLTRTQADPVPRATRLPVEAKDCHIAQFKSLSPKELQAFKKAKDAIEKRLLEKDLRCS  
 SHLFPRAWDLKQLVQERPKALQAEVALTLKVWENMTDSALATILGQPLHTLSHIHSQLQTCTQLQATAE  
 PRSPSRRLSRWLHRLQEAQSKETPGCLEASVTSNLFRLTRDLKCVANGDQCV

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV



[View online »](#)

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9007\\_f12.zip](https://cdn.origene.com/chromatograms/mm9007_f12.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001024673

**ORF Size:** 579 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001024673.2](#), [NP\\_001019844.2](#)

**RefSeq Size:** 582 bp

**RefSeq ORF:** 582 bp

**Locus ID:** 330496


UniProt ID: [Q4VK74](#)

Cytogenetics: 7

MW: 22.2 kDa

**Gene Summary:** Cytokine with antiviral, antitumour and immunomodulatory activities. Plays a critical role in the antiviral host defense, predominantly in the epithelial tissues. Acts as a ligand for the heterodimeric class II cytokine receptor composed of IL10RB and IFNLR1, and receptor engagement leads to the activation of the JAK/STAT signaling pathway resulting in the expression of IFN-stimulated genes (ISG), which mediate the antiviral state. Has a restricted receptor distribution and therefore restricted targets: is primarily active in epithelial cells and this cell type-selective action is because of the epithelial cell-specific expression of its receptor IFNLR1. Seems not to be essential for early virus-activated host defense in vaginal infection, but plays an important role in Toll-like receptor (TLR)-induced antiviral defense. Plays a significant role in the antiviral immune defense in the intestinal epithelium. Exerts an immunomodulatory effect by up-regulating MHC class I antigen expression.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MR224475