

## Product datasheet for **MR224483**

### Ypel4 (NM\_001005342) Mouse Tagged ORF Clone

#### Product data:

Product Type: Expression Plasmids  
 Product Name: Ypel4 (NM\_001005342) Mouse Tagged ORF Clone  
 Tag: Myc-DDK  
 Symbol: Ypel4  
 Mammalian Cell Selection: Neomycin  
 Vector: pCMV6-Entry (PS100001)  
 E. coli Selection: Kanamycin (25 ug/mL)  
 ORF Nucleotide Sequence: >MR224483 ORF sequence  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGCCAGCTGTGACCCTGGCCCGCCCCGCATGTTTACCCACCAAGACTTCCGCAGTTACCTGCCCC  
 GCTGCCACCGCACTTACAGTTGCGTCCACTGTCGAGCACACTTGGCCAAACACGATGAGCTTATTTCCAA  
 GTCCTTCCAAGGAAGCCATGGCCGGGCTACCTGTTTAACTCCGTGGTCAACGTGGTGTGGGCCAGCT  
 GAACAGCGCCTCCTGCTCACAGGACTCCACTCGGTAGCTGACATTTTCTGCGAAAGCTGCAAGACCACAC  
 TGGGCTGGAAATATGAACAGGCTTTTGAGACCAGCCAGAAATACAAAGAAGGGAAGTACATCATTGAGAT  
 GTCACACATGGTAAAGACAACGGCTGGGAT

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR224483 protein sequence  
 Red=Cloning site Green=Tags(s)  
 MPSCDPGPAPACLPTKTFRSYLPRCHRYSVHCRAHLAKHDELISKSFQGSHGRAYLFNSVNVVNGCGPA  
 EQRLLLTGLHSVADIFCESCKTTLGWKYEQAFETSQKYKEGKYIEMSHMVKDNGWD

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM\_001005342

ORF Size: 384 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001005342.2](#), [NP\\_001005342.1](#)

RefSeq Size: 752 bp

RefSeq ORF: 384 bp

Locus ID: 241525  
UniProt ID: [Q65Z93](#)  
Cytogenetics: 2 D  
MW: 14.3 kDa

**Product images:**

Circular map for MR224483