

Product datasheet for **MR224491L4V**

Haus2 (NM_025475) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Haus2 (NM_025475) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Haus2
Synonyms:	1700101G24Rik; 2410007P03Rik; AU020217; Cep27
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_025475
ORF Size:	702 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224491).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025475.3 , NP_079751.1
RefSeq Size:	3170 bp
RefSeq ORF:	705 bp
Locus ID:	66296
UniProt ID:	Q9CQS9
Cytogenetics:	2 E5
Gene Summary:	Contributes to mitotic spindle assembly, maintenance of centrosome integrity and completion of cytokinesis as part of the HAUS augmin-like complex.[UniProtKB/Swiss-Prot Function]



[View online »](#)