

Product datasheet for **MR224509L3V**

Cyp19a1 (NM_007810) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cyp19a1 (NM_007810) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cyp19a1
Synonyms:	Ar; ArKO; Cyp19; Int-5; Int5; p450arom
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_007810
ORF Size:	1509 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224509).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007810.3
RefSeq Size:	2373 bp
RefSeq ORF:	1512 bp
Locus ID:	13075
UniProt ID:	P28649
Cytogenetics:	9 29.49 cM
Gene Summary:	Catalyzes the formation of aromatic C18 estrogens from C19 androgens.[UniProtKB/Swiss-Prot Function]



[View online »](#)