

## Product datasheet for **MR224514**

### **Bloc1s6 (NM\_019788) Mouse Tagged ORF Clone**

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Bloc1s6 (NM_019788) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Bloc1s6
Synonyms:	BLOC-1; pa; Pldn; Stx13bp1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR224514 representing NM_019788 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAGTGTTCCTGAGCCACCGCTCCAGACGGGGTCCTTACGGGGCCTCCGACAGCCTGGAGGCCGGCG  
AGCCACGCCGGGCTTAAGCGACACGTCTCCAGATGAAGGACTGATAGAGGACTTCCTGTGGACGATAG  
AGCTGTAGAGCACTTGGTGGGAGGCCTGCTGTCCCACTACTTACCAGATCTGCAGAGGTCCAAGCGAGCA  
CTCCAGGAACACACAGAACCAAGTTGTGTTACTAGATACACTGGAACAAGAGATTTCAAATTTAAAG  
AATGCCACTCCATGTTGGACATCAATGCTTTGTTACAGAGGCGAAGCACTATCACGCCAAGCTGGTGAC  
TATAAGGAAGGAGATGCTGCTGCCATGAGAAAACGTCCAAGTTGAAAAAAGAGCACTGAAGCTGCAG  
CAGAAGCGGCAGAGAGAAGACTGGAGAGGGAGCAGCAACCGGAAAAGGAATTTGAAAGAGAAAAGCAGC  
TGACAGCCAAGCCAGCCAAGAGGACA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: [https://cdn.origene.com/chromatograms/mm9050\\_c12.zip](https://cdn.origene.com/chromatograms/mm9050_c12.zip)

Restriction Sites: Sgfl-Mlul



[View online »](#)

**Cloning Scheme:**



**ACCN:** NM\_019788

**ORF Size:** 516 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_019788.3](#), [NP\\_062762.1](#)

**RefSeq Size:** 3494 bp

**RefSeq ORF:** 519 bp

**Locus ID:** 18457


**UniProt ID:** [Q9R0C0](#)

**Cytogenetics:** 2 60.66 cM

**MW:** 20.1 kDa

**Gene Summary:** Component of the BLOC-1 complex, a complex that is required for normal biogenesis of lysosome-related organelles (LRO), such as platelet dense granules and melanosomes. In concert with the AP-3 complex, the BLOC-1 complex is required to target membrane protein cargos into vesicles assembled at cell bodies for delivery into neurites and nerve terminals. The BLOC-1 complex, in association with SNARE proteins, is also proposed to be involved in neurite extension. May play a role in intracellular vesicle trafficking, particularly in the vesicle-docking and fusion process.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MR224514