

Product datasheet for **MR224544L3V**

Dcstamp (NM_029422) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dcstamp (NM_029422) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dcstamp
Synonyms:	4833414I07Rik; DC-STAMP; FIND; mDC-STAMP; Tm7sf4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_029422
ORF Size:	1410 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224544).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_029422.3 , NP_083698.1
RefSeq Size:	1993 bp
RefSeq ORF:	1413 bp
Locus ID:	75766
UniProt ID:	Q7TNJ0
Cytogenetics:	15 B3.1



[View online »](#)

Gene Summary:

Probable cell surface receptor that plays several roles in cellular fusion, cell differentiation, bone and immune homeostasis. Plays a role in TNFSF11-mediated osteoclastogenesis. Cooperates with OCSTAMP in modulating cell-cell fusion in both osteoclasts and foreign body giant cells (FBGCs). Participates in osteoclast bone resorption. Involved in inducing the expression of tartrate-resistant acid phosphatase in osteoclast precursors. Plays a role in haematopoietic stem cell differentiation of bone marrow cells toward the myeloid lineage. Inhibits the development of neutrophilic granulocytes. Plays also a role in the regulation of dendritic cell (DC) antigen presentation activity by controlling phagocytic activity. Involved in the maintenance of immune self-tolerance and avoidance of autoimmune reactions.
[UniProtKB/Swiss-Prot Function]