

Product datasheet for **MR224576L4V**

Josd2 (NM_025368) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Josd2 (NM_025368) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Josd2
Synonyms:	1110007C05Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_025368
ORF Size:	564 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224576).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025368.4 , NP_079644.1
RefSeq Size:	853 bp
RefSeq ORF:	567 bp
Locus ID:	66124
UniProt ID:	Q9CR30
Cytogenetics:	7 B3
Gene Summary:	Cleaves 'Lys-63'-linked poly-ubiquitin chains, and with lesser efficiency 'Lys-48'-linked poly-ubiquitin chains (in vitro). May act as a deubiquitinating enzyme (By similarity). [UniProtKB/Swiss-Prot Function]


[View online »](#)