

Product datasheet for **MR224586A1V**

Mouse Fxn (NM_008044) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Fxn (NM_008044) AAV Particle
Tag: Myc-DDK
Symbol: Fxn
Synonyms: FA; FARR; Frda; X25
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR224586 representing NM_008044
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTGGGCGTTTCGGAGGTCGCGCAGCCGTGGGCTTGCTGCCCCGACGGCGTCCCGGGCCTCCGCCTGGG
TCGGGAACCCGCGCTGGAGGGAACCGATCGTAACCTGCGGCCGCCGAGGCCTACATGTCACAGTCAACGC
CGGCGCCACCCGCCACGCCCATTTGAACCTCCACTACCTCCAGATTCTGAACATCAAAAAGCAGAGCGTC
TGCGTGGTGCATTTGAGGAACCTGGGGACATTGGACAACCAAGCTCTCTAGACGAGACAGCGTATGAAA
GACTGGCGGAAGAGACCCTGGACTCCCTGGCCGAGTTCTTTGAAGACCTCGCAGACAAGCCCTATACCTT
GGAGGACTACGATGTCTCTTTTGGGGATGGCGTGCTCACCATTAAAGCTGGGCGGGGATCTAGGGACCTAC
GTGATCAACAAGCAGACCCCAACAAGCAAATCTGGCTGTCTTCTCCTCCAGCGGCCCAAGCGCTATG
ACTGGACCGGGAAGAACTGGGTGTAATCTCATGACGGCGTGCTCTGCATGAGCTGCTGGCCAGGGAGCT
GACTAAAGCTTTAAACACCAAACTGGACTTGTCTTCATTGGCCTATTCTGGAAAAGGCACT

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR224586 representing NM_008044
Red=Cloning site Green=Tags(s)

MWAFGGRAAVGLLPRTASRASAWVGNPRWREPIVTCGRRGLHVTNAGATRHAHLNLHYLQILNIKKQSV
CVVHLRNLGTLNPPSLDETAYERLAEETLDSLAEFFEDLADKPYTLEDYDVVSFGDGLTIKLGDLGT
VINKQTPNKQIWLSSPSSGPKRYDWTGKNWVYSHDGVSLHELLARELTKALNTKLDLSSLAYSGKGT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_008044
ORF Size:	621 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_008044.2</u>
RefSeq Size:	1095 bp
RefSeq ORF:	624 bp
Locus ID:	14297
UniProt ID:	<u>O35943</u>
Cytogenetics:	19 B
MW:	23.4 kDa