

Product datasheet for **MR224615L3V**

Kcnj13 (NM_001110227) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Kcnj13 (NM_001110227) Mouse Tagged ORF Clone Lentiviral Particle |
| Symbol: | Kcnj13 |
| Synonyms: | 100040591 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_001110227 |
| ORF Size: | 1080 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR224615). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001110227.2 , NP_001103697.1 |
| RefSeq Size: | 3056 bp |
| RefSeq ORF: | 1083 bp |
| Locus ID: | 100040591 |
| UniProt ID: | P86046 |
| Cytogenetics: | 1 |



[View online »](#)

Gene Summary:

Inward rectifier potassium channels are characterized by a greater tendency to allow potassium to flow into the cell rather than out of it. Their voltage dependence is regulated by the concentration of extracellular potassium; as external potassium is raised, the voltage range of the channel opening shifts to more positive voltages. The inward rectification is mainly due to the blockage of outward current by internal magnesium. KCNJ13 has a very low single channel conductance, low sensitivity to block by external barium and cesium, and no dependence of its inward rectification properties on the internal blocking particle magnesium (By similarity).[UniProtKB/Swiss-Prot Function]