

## Product datasheet for **MR224637L3V**

### Lmtk2 (NM\_001081109) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Lmtk2 (NM_001081109) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lmtk2
Synonyms:	2900041G10Rik; A330101P12Rik; aatyk2; BREK; cprk; KPI-2; KPI2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001081109
ORF Size:	4413 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224637).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001081109.1</a>
RefSeq Size:	8114 bp
RefSeq ORF:	4416 bp
Locus ID:	231876
UniProt ID:	<a href="#">Q3TYD6</a>
Cytogenetics:	5 G2
Gene Summary:	Phosphorylates PPP1C, phosphorylase b and CFTR.[UniProtKB/Swiss-Prot Function]


[View online »](#)