

Product datasheet for **MR224668L3V**

Ei24 (NM_007915) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ei24 (NM_007915) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ei24
Synonyms:	AA536736; AI115355; PIG8
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_007915
ORF Size:	1074 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224668).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007915.5 , NP_031941.1
RefSeq Size:	2108 bp
RefSeq ORF:	1077 bp
Locus ID:	13663
UniProt ID:	Q61070
Cytogenetics:	9 20.68 cM
Gene Summary:	Acts as a negative growth regulator via p53-mediated apoptosis pathway. Regulates formation of degradative autolysosomes during autophagy (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)