

Product datasheet for **MR224696**

Trmt1 (NM_001164560) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Trmt1 (NM_001164560) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Trmt1
Synonyms:	6720406L13Rik; 6720477L24; D8Ertd812e
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:


>MR224696 representing NM_001164560
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCGAGAAGCATGCACCGTGA CTCTGGATCCTCTGTAGAGCTCAGTTTCATGGAACGGAAGGCC
 AGAAGCCTCCCTCTCCACCAGCGATGGAGAACGGCACCCGGTCTGGGAAGAGCGCCACCCGCGGACCC
 GGTTCGACTGTCACTGAGGGCGCCGCAAAATCGTCTTTCCAGTGCTAACGAGTTTTCTATAACCCG
 GTGCAGGAATCAACCGGACCTGACGTGTCTGTGATCACAGAGTTTGTCTGAATTCATCTTGGGGCCA
 AGGGAATCCAGATCAAGTTCCAGGAGAGAAGGACTCAGAGAAAAAGCTGTGGACTTGTCTGATCAAGA
 AGAGGAAACAGCTGAAAGAATGAAATCTGGCCCCAGGGGACTGGCCACGGACTGCTGCAGTGGGTGAG
 ATCTGTGAGGAAGGTCTACGAGTCTGGAAGGCCTGGCAGCTTCCGGTCTGCGCTCCATCCGTTTTGCC
 TGGAGTTCTGGCCTCCAGTCTGTGGTGGCAATGATGCCTCTGCCCGGCTGTAGAGCTTATGCACAG
 AAATGTGGAGCTCAATGGTGTGGCAGATCTTGTTCAGCCCAACCAGGCAGATGCCCGGATGCTGATGAC
 CAACACAAAAGGCACCTGAGCGATTTGATGTCATTGACCTGGACCCCTATGGCAGCCCCGCCCTTCC
 TGGATGCAGCAGTGCAGGCTGTGAGTGTGGAGGACTTCTGTGTGTACCTGCACAGACATGGCAGTCTCT
 GGCGGGAAACAGCGGAGAGACATGCTACAGCAAGTATGGGGCCATGGCCCTCAAGAGCCGAGCCTGCCAT
 GAAATGGCCCTGCGGATTGTGCTGCACAGCCTGGACCTGCACGCCAAGTCTACCAGCGCTATATCGTGC
 CGCTGCTCAGCATCAGTGTGACTTCTATATCCGGTCTTTGTACGTGTTTTACAGGCCAGGCCAAGGT
 TAAGTCTCAGCCAGCAAACAGGCCCTGGTGTCCAGTGTGGGTTGTGGAGCCTTCTACCTGCAACGC
 CTTGGCAAAGCATCAGGAGACCTGGTGGCAGGATCAAGTCTCTGCAGCCTGTGGTCTCCAGTAACCT
 CCGAGTGTGAGCACTGTGGACAGAGACACCAGCTTGGGGCCCCATGTGGGCAGAGCCCATCCATGACCT
 GGACTTTGTGGCCGAGTTCTTGACGCAGTGACCACGAATCCTGGCCGCTTCCACACTTCCATGAGGATC
 CAAGGCGTCTCAGTGTGTTACTGAGGAACTTCTGATGTCCTCTCTACTACACACTGGACCAACTGA
 GCAGCACTATCCACTGCAACACACCCCGCTTCTGCAGCTTCGGTCTCAGCTCTCCTTCATGCTGGTTCCG
 AGTCTCTCTCTCCATGCTTGTAAAGATGCAGTGAAGACAGATGCTCCCCCTGAGGCCCTCTGGGACATC
 ATGCGTTGCTGGGAAAAGGAATGTCCAGTGAACGGGAGCGGCTGTCTGAAAGCAGTCCAGCATTCCGAA
 TTCTTGTGTGGAGCCAGGCTGAAGGCCAAGTCAACATCCGGGAAGACGCCAACCCAGCTCCCGCCA
 GCGAGGACTTAAGCGCTTCCAGGCCAATCCAGAAGCCAAGTGGGGTCCCCGGCCCGTCCCGGCCAGGG
 GGCAAGGCAGCAAGTGAAGATCTAGCAGGGAGGCGGAGGCTGCTACAGAATAAGAGGAAGGAGCCAGCTG
 AGGACCCGGCCAGCGTGCAGCTCGCCTTAAGACATTTCTTGCAAACGATTCAAAGAGGGTACCTGTCA
 ACTAGGGGACCAGTGTGTTATTACACAGTCTGCGAGCTCTGTGGCCTCTGGTGTATATCCCCATCGAA
 GAATGTCCAGAGACTACCACCAAGATCTCCCCTGGACCAAGGCTGCTGTGGGGTATTCTGGGCCAG
 GTGTAGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Plasmid Map:



ACCN: NM_001164560

ORF Size: 1968 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001164560.1](#), [NP_001158032.1](#)

RefSeq Size: 2083 bp

RefSeq ORF: 1971 bp

Locus ID: 212528

UniProt ID: [Q3TX08](#)

Cytogenetics: 8 41.02 cM

MW: 72.1 kDa

Gene Summary: Dimethylates a single guanine residue at position 26 of most tRNAs using S-adenosyl-L-methionine as donor of the methyl groups.[UniProtKB/Swiss-Prot Function]