

Product datasheet for **MR224741**

Zbtb44 (NM_001115130) Mouse Tagged ORF Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | Zbtb44 (NM_001115130) Mouse Tagged ORF Clone |
| Tag: | Myc-DDK |
| Symbol: | Zbtb44 |
| Synonyms: | 6030404E16Rik; Btbd15; E230011F09 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |



[View online »](#)

ORF Nucleotide Sequence:

>MR224741 representing NM_001115130
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGGTGTGAAAACATTTACTCACAGCTCCTCTCCACAGCCAGGAAATGCTGGGAAACTAAATATGC
 TCCGAAATGATGGACATTTTTGTGATATCACTATTCGTGTCCAGGACAAAATCTCCGGCTCACAGGT
 AGTACTAGCAGCTTGCAGTGATTCTTTTCGAACCAAATTTGTTGGCCAAACAGAAGATGAGAACAAAAAT
 GTGTTGGATCTGCATCATGTTACAGTACTGGCTTTATACCTCTTTTAGAATACGCTTACACAGCCACTC
 TGTCAATTAACACGAAAATATTATTGATGTTCTAGCGGCAGCCAGCTATATGCAAATGTTTAGTGTTC
 TAGTACATGCTCAGAGTTCATGAAGTCAAGTATTTTGTGGAATACACCTAATAGCCAACCCAGAAAAGAGT
 CTAGATGCTGGACAAGAAAATAGCTCCAAGTCAATTTACTTCTCGAGATGGAAGCATTCTCCAGTGT
 CTTAGAGTGCAGTGTGTGGAGAGGACCATTCTGTCTGCCGAGAGTCCCGGAGGAAGCGCAAAGCTA
 CATTGTTATGTCTCCCGAGAGTCCCGTGAATGTAGCACCAAACAAGTCCCCCGAGTACTGAATTCT
 TCAGCTTCTATGCAGAAAATAGAAGTCAGCCAGTTGACTCTTCTTTAGCTTTTCCCTGGACTTTCCCTT
 TTGGAATAGATCGAAGGATTCAGCCTGAGAAAGCAAAGCAGGCAGAAAATACCCGGACTTTAGAATTGCC
 TGGCCATCCGAGGCAGGTAGAAGAGTGGCTGACTATGTGACTTGTGAGAGCACAAAACCCACTTGCTT
 CTAGGTACAGAAGAAGATGTCCGGGTCAAAGTAGAAAAGATTAAGTGATGAGGAGGTCATGAGGAAGTGT
 CCCAGCCTGTTAGTGCATCTCAGAGTCACTGAGTATCAGCAGACAGTCCCGGAAGTGAACAGTTCA
 AGAGGACCTTCTGATCAGTCCACAGTCTTCTCTATAGGCTCAGTAGAAGGAGTCACTGAAGGATTG
 CCTACACTCAGAGCACCTTAGCACCAATGCTCATGCAGATGACGATGACAGATTGGAAAATGTTCAAGT
 ATCCCTATCAACTCTACATTGCTCTTCTACCAGCAGTACAGAGCGGCAAGTCCAAATGGTCCAGACAG
 ACCGTTTTAGTGTCCAACCTGTGGGGTCCGCTTACCCGTATTCAGAATCTAAAACAGCACATGCTTATC
 CACTCAGGAATTAACCATTTTCAGTGTGACTGCTGTGAAAAAAGTTTACCAGGGCTTACTCGCTAAAGA
 TGCATCGCCTAAAGCACGAAGTATTTCT

ACGCGTACGCGGCCGCTCGAGCAGAAAATCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR224741 representing NM_001115130
 Red=Cloning site Green=Tags(s)

MGVKFTFHSSSSHSQEMLGKLNMLRNDGHFCDITIRVQDKIFRAHKVVLAAACSDFFRTRKLVGQTEENKN
 VLDLHHVTVTGFIPLLEYAYTATLSINTENIIDVLAASYSYQMFVASTCSEFMKSSILWNPNSQPEKS
 LDAGQENSSNCFSTRDGSISPVSECSAERTIPVCRESRRKRKSYIVMSPEPVCSTQTSSPQVLNS
 SASYAENRSQPVDSSLAFPWTFPFIDRRIQPEKAKQAENTRTLELPGPSEAGRRVADYVTCSTKPTLP
 LGTEEDVRVKVERLSDEEVHEEVSPVSAQSSSDQQTVPGSEPVQEDLLISQSSSIGSVDEGVTEGL
 PTLQSTSTNAHADDDRLNENVQYPYQLYIAPSTSSSTERPSNPGPDRPFQCPCTCGVRFTRIQNLKQHMLI
 HSGIKPFQDCCKKFFTRAYSCLKMHLKHEVIS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/ja2394_b12.zip

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: NM_001115130

ORF Size: 1362 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)


OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001115130.1](#), [NP_001108602.1](#)
 RefSeq Size: 8882 bp
 RefSeq ORF: 1362 bp
 Locus ID: 235132
 UniProt ID: [Q8R0A2](#)
 Cytogenetics: 9 A4
 MW: 50.2 kDa

Product images:



Circular map for MR224741