

Product datasheet for **MR224747L4V**

Eif5a2 (NM_177586) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Eif5a2 (NM_177586) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Eif5a2
Synonyms:	9630038B20; eIF5AII
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_177586
ORF Size:	459 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224747).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_177586.5
RefSeq Size:	4992 bp
RefSeq ORF:	462 bp
Locus ID:	208691
UniProt ID:	Q8BGY2
Cytogenetics:	3 A3



[View online »](#)

Gene Summary:

mRNA-binding protein involved in translation elongation. Has an important function at the level of mRNA turnover, probably acting downstream of decapping. Involved in actin dynamics and cell cycle progression, mRNA decay and probably in a pathway involved in stress response and maintenance of cell wall integrity. Functions as a regulator of apoptosis. Mediates effects of polyamines on neuronal process extension and survival. May play an important role in brain development and function, and in skeletal muscle stem cell differentiation (By similarity).[UniProtKB/Swiss-Prot Function]