

Product datasheet for **MR224756**

Il17re (NM_145826) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Il17re (NM_145826) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Il17re
Synonyms:	AA589509; IL-17RE; Il25r
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR224756 representing NM_145826
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGGGAGCCCCAGACTGGCAGCCTTGCTCCTGTCTCTCCGCTACTGCTCATCGGCCTCGCTGTGTCTG
 CTGGGTTGCCTGCCCTGCCTGCGGAGTTGGACCAGCCACTGTCTCCTGGCCTACCGTGTGGATAAACG
 TTTTGCTGGCCTTCAGTGGGGCTGGTCCCTCTCTTGGTGGAAATCTAAAAGTCTCCTAAATTTGAA
 GACTATTGGAGGCACAGGACACCAGCATCCTTCCAGAGGAAGCTGCTAGGCAGCCCTTCCCTGTCTGAGG
 AAAGCCATCGAATTTCCATCCCCCTCAGCCATCTCCACAGAGGCCAACGCACCAAAAGGGCCAGCC
 TTCAGTGCAGAAGGAAGAGAACATCTCCCTGAAGCAGGGTCACAAAAGTGTGGAGGACCTGAATTTCC
 TTTGATTTGCTGCCCGAGGTGCAGGCTGTTCCGGTGACTATTCTGCAGGCCCAAGGCCAGTGTGCGCC
 TTTGTTATCAGTGGGCACTGGAATGTGAAGACTTGAGTAGCCCTTTTGATACCCAGAAAATTGTGTCTGG
 AGGCCACACTGTAGACCTGCCTTATGAATTCCTTCTGCCCTGCATGTGCATAGAGGCCTCCTACCTGCAA
 GAGGACACTGTGAGGCGCAAAAAGTGTCCCTTCCAGAGCTGGCCTGAAGCTTATGGCTCAGACTTCTGGC
 AGTCAATACGCTTCACTGACTACAGCCAGCACAATCAGATGGTCAATGGCTCTGACTCCGCTGCCCACT
 GAAACTGGAGGCCTCCCTCTGCTGGAGGCAGGCCACTCACACCCTGCGAAACCCTTCCCAACGCCACA
 GCACAGGAGTCAGAAGGATGGTATATCCTGGAGAATGTGGACTTGACCCCCAGCTCTGCTTTAAGTTCT
 CATTTGAAAACAGCAGCCACGTTGAATGTCCCCACCAGAGTGGCTCTCTCCATCCTGGACTGTGAGCAT
 GGATACCCAGGCCAGCAGCTGACGCTTCACTTTTCTTCGAGGACATATGCCACCTTCAGTGTGCTGCTGG
 AGTGAACAGGTTGGGGCCGGATACCCCATGCCTCCTGTGTACAGCATCAGCCAGACCCAGGGCTCAG
 TCCCAGTACGCTAGACCTCATCATCCCCTTCCCTGAGGCAGGAGAATTGCATCCTGGTGTGGAGGTCAGA
 TGTCCATTTTGCCTGGAAGCACGTCTTGTGTCTGATGTCTCCATAGACACCTCGGGCTCTTGATCCTG
 GCACTGCTGGCTCTACCGCTCTAGTGGGTGTAGTTCTGGTCTCCTCGGCCGGCGCCTACTGCCAGGCT
 CCGGTCGAACAAGGCCAGTTTTACTCTACATGCAGCGGACTCAGAGGCACAGCGACGCCTGGTGGGAGC
 TTTGGCCGAAGTCTGCGGACGGCGCTGGGAGGTGGACGCGACGTGATCGTGGATCTCTGGGAAGGGACG
 CACGTAGCACGCATTGGACCACTGCCGTGGCTTTGGGCAGCGCGGGAGCGCGTGGCGCGGGAGCAGGGCA
 CAGTGTGCTCCTGTGGAAGTGTGCGGGTCCCAGCACCGCCTGCAGCGGTGACCCGAGGCTGCGTCCCT
 TCGCACCTTGTGTGCGCTGCTCCACGTCCGCTGCTGCTCGCCTACTTCAGTCGCCTCTGCGCAAAGGT
 GACATCCCCCGGCCGCTGCGCGCTGCCACGCTACCGCCTGCTTCGTGACCTGCCGCGCCTGCTGAGAG
 CACTGGATGCTCAGCCTGCCACCCTAGCCTCCAGCTGGAGTACCTTGGGGCTAAGCGGTGCTTGAAAAA
 CCGTCTGGAGCAGTGTACCTGTGGAAGTGGAGGCTGCCAAAGATGACTACCAAGGCTCAACCAATAGT
 CCTGTGGTTTCAGCTGTCTG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR224756 representing NM_145826
Red=Cloning site Green=Tags(s)


MGSPRLAALLLSLPLLLIGLAVSARVACPCLRSWTSHELLAYRVDKRFAGLQWGWFPLLVRKSKSPPKFE
 DYWRHRTPASFQRKLLGSPSLSEESHRISSIPSSAISHRGQRTKRAQPSAAEGREHLPEAGSQKCGGPEFS
 FDLLPEVQAVRVTIPAGPKASVRLCYQWALECEDLSSPFDTKIIVSGGHTVDLPYEFLLPCMIEASYLQ
 EDTVRRKKCPFQSWPEAYGSDFWQSIRFTDYSQHNQVMALTLRCPLKLEASLCWRQDPLTPCETLPNAT
 AQESEGWYILENVDLHPQLCFKFSFENS SHVECPHQSGSLPSWTVSMDTQAQQLTLHFSSRTYATFSAAW
 SDPGLGPDTPMPVYSISQTQGSVPVTLDLIIPFLRQENCILVWRSDVHFAWKHVLCPDVSHRHLGLLIL
 ALLALTALVGVVLLGRLLPGSGRTRPVLLHAADSEAQRRLVGALAE LLRTALGGGRDIVDLWEGT
 HVARIGPLPWLWAARERVAREQGTVLLWNCAGPSTACSGDPQAASLRTLLCAAPRPLLLAYFSRLCAKG
 DIPRPLRALPRYRLRLDLPRLLRALDAQPATLASSWSHLGAKRCLKNRLEQCHLLEEAAKDDYQGSTNS
 PCGFSC L

TRTRPLEQKLISEEDLANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9005_a11.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_145826

ORF Size: 1911 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_145826.6](#)

RefSeq Size: 2933 bp

RefSeq ORF: 1914 bp

Locus ID: 57890


UniProt ID: [Q8BH06](#)

Cytogenetics: 6 E3

MW: 71.3 kDa

Gene Summary: Specific functional receptor for IL17C, signaling through the NF-kappa-B and MAPK pathways. Requires TRAF3IP2 /ACT1 for signaling. Crucial regulator in innate immunity to bacterial pathogens, such as *Citrobacter rodentium*. Isoform 4 and isoform 5 may be either cytoplasmic inactive or dominant active forms. Isoform 2 and isoform 3 may act as soluble decoy receptors.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR224756