

# Product datasheet for MR224787L2

### Tnfsf12 (NM\_011614) Mouse Tagged Lenti ORF Clone

#### **Product data:**

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

| Product Type:                | Expression Plasmids   |
|------------------------------|---|
| Product Name:                | Tnfsf12 (NM_011614) Mouse Tagged Lenti ORF Clone  |
| Tag:                         | mGFP  |
| Symbol:                      | Tnfsf12   |
| Synonyms:                    | Apo3l; Dr3l; Dr3lg; Tweak   |
| Mammalian Cell<br>Selection: | None  |
| Vector:                      | pLenti-C-mGFP (PS100071)  |
| E. coli Selection:           | Chloramphenicol (34 ug/mL)  |
| ORF Nucleotide<br>Sequence:  | The ORF insert of this clone is exactly the same as(MR224787).                            |
| <b>Restriction Sites:</b>    | Sgfl-Mlul   |
| Cloning Scheme:              |   |
|                              | Cloning sites used for ORF Shuttling:<br>Sgf I ORF Mlu I<br>GCG ATC GCC ATG// NNN ACG CGT |



ACCN: ORF Size: NM\_011614 747 bp




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US


|                     | sf12 (NM_011614) Mouse Tagged Lenti ORF Clone – MR224787L2   |
|---------------------|--|
| OTI Disclaimer:     | Due to the inherent nature of this plasmid, standard methods to replicate additional amounts<br>of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore,<br>OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts<br>of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a<br>reduced cost. Please contact our customer care team at <u>custsupport@origene.com</u> or by<br>calling 301.340.3188 option 3 for pricing and delivery. |
|                     | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>  |
| OTI Annotation:     | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| Components:         | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).   |
| Reconstitution Meth | <ul> <li>centrifuge at 5,000xg for 5min.</li> <li>Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>Close the tube and incubate for 10 minutes at room temperature.</li> <li>Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ul>   |
| RefSeq:             | <u>NM 011614.2</u>   |
| RefSeq Size:        | 1515 bp  |
| RefSeq ORF:         | 750 bp   |
| Locus ID:           | 21944  |
| UniProt ID:         | <u>054907</u>  |
| Cytogenetics:       | 11 B3  |
| Gene Summary:       | Binds to FN14 and possibly also to TNRFSF12/APO3. Weak inducer of apoptosis in some cell types. Mediates NF-kappa-B activation. Promotes angiogenesis and the proliferation of endothelial cells. Also involved in induction of inflammatory cytokines. Promotes IL8 secretion (By similarity).[UniProtKB/Swiss-Prot Function]   |

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

## **Product images:**



Circular map for MR224787L2



Double digestion of MR224787L2 using Sgfl and Mlul

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US