

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Product datasheet for MR224789L4V

## Zdhhc17 (NM\_172554) Mouse Tagged ORF Clone Lentiviral Particle

## **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Zdhhc17 (NM_172554) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Zdhhc17
Synonyms:	A230053P19Rik; BB187739; D130071N24Rik; Hip14
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_172554
ORF Size:	1896 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224789).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<u>NM 172554.2, NP 766142.2</u>
RefSeq Size:	4563 bp
RefSeq ORF:	1899 bp
Locus ID:	320150
UniProt ID:	<u>Q80TN5</u>
Cytogenetics:	10 D1



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	Zdhhc17 (NM_172554) Mouse Tagged ORF Clone Lentiviral Particle – MR224789L4V	
Gene Summarv:	Palmitovltransferase specific for a subset of neuronal proteins, including SNAP25.	

Palmitoyltransferase specific for a subset of neuronal proteins, including SNAP25, DLG4/PSD95, GAD2, SYT1 and HTT (PubMed:15489887, PubMed:15603741, PubMed:25253725). May be involved in the sorting or targeting of critical proteins involved in the initiating events of endocytosis at the plasma membrane (By similarity). May play a role in Mg(2+) transport (By similarity).[UniProtKB/Swiss-Prot Function]

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US