

Product datasheet for **MR224795L3V**

PPP3R2 (NM_001004025) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ppp3r2 (NM_001004025) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ppp3r2
Synonyms:	CaNB2; CnB2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001004025
ORF Size:	540 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224795).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001004025.4 , NP_001004025.1
RefSeq Size:	2985 bp
RefSeq ORF:	540 bp
Locus ID:	19059
UniProt ID:	Q63811
Cytogenetics:	4 B1
Gene Summary:	Regulatory subunit of calcineurin, a calcium-dependent, calmodulin stimulated protein phosphatase. Confers calcium sensitivity.[UniProtKB/Swiss-Prot Function]



[View online »](#)