

Product datasheet for **MR224796L4V**

Tnnt2 (NM_001130181) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tnnt2 (NM_001130181) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tnnt2
Synonyms:	cTnT; Tnt
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001130181
ORF Size:	957 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224796).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001130181.2 , NP_001123653.4
RefSeq Size:	1165 bp
RefSeq ORF:	960 bp
Locus ID:	21956
UniProt ID:	P50752
Cytogenetics:	1 59.32 cM
Gene Summary:	Troponin T is the tropomyosin-binding subunit of troponin, the thin filament regulatory complex which confers calcium-sensitivity to striated muscle actomyosin ATPase activity. [UniProtKB/Swiss-Prot Function]



[View online »](#)