

Product datasheet for **MR224814L4V**

Dido1 (NM_011805) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dido1 (NM_011805) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dido1
Synonyms:	6720461J16Rik; C130092D22Rik; D130048F08Rik; Datf; DATF-1; Datf1; di; dido; DIO; DIO-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_011805
ORF Size:	1842 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224814).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011805.2 , NP_035935.2
RefSeq Size:	4662 bp
RefSeq ORF:	1845 bp
Locus ID:	23856
UniProt ID:	Q8C9B9
Cytogenetics:	2 H4



[View online »](#)

Gene Summary:

This gene encodes a transcription factor involved in apoptosis. The encoded protein functions in cell cycle progression and plays a role in chromosomal stability. This protein regulates the self-renewal of embryonic stem cells. Disruption of this gene in mice causes symptoms similar to myelodysplastic/myeloproliferative diseases in humans. Mice lacking this gene show severely reduced fertility. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2014]