

Product datasheet for **MR224895L4V**

Scgb3a1 (NM_054037) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Scgb3a1 (NM_054037) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Scgb3a1
Synonyms:	HIN; HIN-1; LuLe; LuLeu2; Pnsp2; UGR; UGRP2; Utgrp2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_054037
ORF Size:	315 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224895).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_054037.1 , NP_473378.1
RefSeq Size:	540 bp
RefSeq ORF:	315 bp
Locus ID:	68662
UniProt ID:	Q920D7
Cytogenetics:	11 29.75 cM



[View online »](#)

Gene Summary:

This gene encodes a small secreted protein that may function in the inhibition of cell growth. Expression of the related gene in humans is downregulated via DNA methylation in breast cancers and other tumor types. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jan 2013]