

Product datasheet for **MR224909L4V**

Brd4 (NM_020508) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Brd4 (NM_020508) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Brd4
Synonyms:	Brd5; HUNK1; MCAP; WI-11513
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_020508
ORF Size:	4200 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224909).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_020508.3 , NP_065254.3
RefSeq Size:	5948 bp
RefSeq ORF:	4203 bp
Locus ID:	57261
UniProt ID:	Q9ESU6
Cytogenetics:	17 17.39 cM
Gene Summary:	This gene was temporarily named bromodomain-containing 5 (Brd5) and was renamed bromodomain-containing 4 (Brd4). [provided by RefSeq, Jul 2008]



[View online »](#)