

## Product datasheet for **MR224915L3V**

### **Ablim2 (NM\_001177697) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Ablim2 (NM_001177697) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ablim2
Synonyms:	A1606905; C230091L11
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001177697
ORF Size:	1851 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224915).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001177697.1</a> , <a href="#">NP_001171168.1</a>
RefSeq Size:	3585 bp
RefSeq ORF:	1854 bp
Locus ID:	231148
UniProt ID:	<a href="#">Q8BL65</a>
Cytogenetics:	5 B3
Gene Summary:	May act as scaffold protein. May stimulate ABRA activity and ABRA-dependent SRF transcriptional activity.[UniProtKB/Swiss-Prot Function]



[View online »](#)