

Product datasheet for **MR224931L3V**

Elf2 (NM_023502) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Elf2 (NM_023502) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Elf2
Synonyms:	2610036A20Rik; A230104O07Rik; AW111824; BB183398; EU32; NERF; NERF-1A; NERF-1B; NERF-2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_023502
ORF Size:	1779 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224931).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_023502.1 , NP_075991.1
RefSeq Size:	2688 bp
RefSeq ORF:	1782 bp
Locus ID:	69257
UniProt ID:	Q9JHC9
Cytogenetics:	3 C
Gene Summary:	Probably transcriptionally activates the LYN and BLK promoters and acts synergistically with RUNX1 to transactivate the BLK promoter.[UniProtKB/Swiss-Prot Function]



[View online »](#)