

## Product datasheet for **MR224956L3V**

### Espn (NM\_207689) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Espn (NM_207689) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Espn
Synonyms:	je
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_207689
ORF Size:	1356 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224956).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_207689.2</a> , <a href="#">NP_997572.1</a>
RefSeq Size:	1647 bp
RefSeq ORF:	1359 bp
Locus ID:	56226
UniProt ID:	<a href="#">Q9ET47</a>
Cytogenetics:	4 82.9 cM



[View online »](#)

**Gene Summary:**

Multifunctional actin-bundling protein. Plays a major role in regulating the organization, dimension, dynamics and signaling capacities of the actin filament-rich microvilli in the mechanosensory and chemosensory cells (PubMed:14657236, PubMed:15190118). Required for the assembly and stabilization of the stereociliary parallel actin bundles. Plays a crucial role in the formation and maintenance of inner ear hair cell stereocilia (PubMed:21455486). Involved in the elongation of actin in stereocilia (PubMed:19287378, PubMed:22264607). In extrastricular hair cells, required for targeting MYO3B to stereocilia tips, and for regulation of stereocilia diameter and staircase formation (PubMed:26926603).[UniProtKB/Swiss-Prot Function]