

## Product datasheet for **MR224979L4V**

### Lrp12 (NM\_172814) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Lrp12 (NM_172814) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lrp12
Synonyms:	A1848829; C820005L12Rik; LRP-12
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_172814
ORF Size:	2574 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224979).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_172814.3</a>
RefSeq Size:	4064 bp
RefSeq ORF:	2577 bp
Locus ID:	239393
UniProt ID:	<a href="#">Q8BUJ9</a>
Cytogenetics:	15 B3.1
Gene Summary:	Probable receptor, which may be involved in the internalization of lipophilic molecules and/or signal transduction. May act as a tumor suppressor (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)