

## Product datasheet for MR224999L3

### Adamts13 (NM\_001001322) Mouse Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Adamts13 (NM_001001322) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Adamts13
Synonyms:	ADAM-TS13; ADAMTS-13; Gm710; vWF-CP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR224999).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.


ACCN:	NM_001001322
ORF Size:	4278 bp



[View online »](#)

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">NM_001001322.2</a> , <a href="#">NP_001001322.1</a>
<b>RefSeq Size:</b>	4580 bp
<b>RefSeq ORF:</b>	4281 bp
<b>Locus ID:</b>	279028
<b>UniProt ID:</b>	<a href="#">Q769J6</a>
<b>Cytogenetics:</b>	2 A3
<b>Gene Summary:</b>	This gene encodes a member of "a disintegrin and metalloproteinase with thrombospondin motifs" (ADAMTS) family of multi-domain matrix-associated metalloendopeptidases that have diverse roles in tissue morphogenesis and pathophysiological remodeling, in inflammation and in vascular biology. In certain mouse strains (C57BL/6, for example) an intracisternal A-type particle (IAP) retrotransposon sequence is located in the intron 23 that causes an alternate splicing event resulting in a shorter transcript variants encoding shorter isoforms. The encoded preproprotein undergoes proteolytic processing to generate an active enzyme that cleaves von Willebrand factor (VWF) in circulating blood. [provided by RefSeq, Jul 2016]

## Product images:



Circular map for MR224999L3