

## Product datasheet for **MR225007L3V**

### Hao1 (NM\_010403) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Hao1 (NM_010403) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hao1
Synonyms:	GOX; Gox1; Hao-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_010403
ORF Size:	1113 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225007).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_010403.2</a> , <a href="#">NP_034533.1</a>
RefSeq Size:	2029 bp
RefSeq ORF:	1113 bp
Locus ID:	15112
UniProt ID:	<a href="#">Q9WU19</a>
Cytogenetics:	2 65.65 cM
Gene Summary:	Has 2-hydroxyacid oxidase activity. Most active on the 2-carbon substrate glycolate, but is also active on 2-hydroxy fatty acids, with high activity towards 2-hydroxy palmitate and 2-hydroxy octanoate (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)