

Product datasheet for **MR225011A1V**

Mouse Zfp784 (NM_001039532) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Zfp784 (NM_001039532) AAV Particle
Tag:	Myc-DDK
Symbol:	Zfp784
Synonyms:	C730049F21Rik; Znf784
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR225011 representing NM_001039532 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCCGCCGCGCGCCCGGATCCCCGATTCCGAGCTACCGACCCGGGAGTCGCCATCCCCGGAGCCGC
CGGACCTGGTCCTGGTGCCTGATGGTCGCCCTGTCACACCCAGGTGACCTCATCGAGATCCAGGTGGT
AAAAGTCACAGATACCACATCGGTCCCTGAGCCCCAGAGCCAGGATCTTTTCATTGTGCCCTGTGTCCT
GCTGCCCTCCGTCTAGTCTCAGAGCTGCTGTTTCATGAGCATGGCCACCTGGCAAGCATGGAGGGCTTG
GGCAGGATGGGATCCTAGCCGTTGTCATGTGTGTGGCCACAGCTGTCCAGGTCCTGCTAGCCTGCGTGC
CCATTACAGCTTGCACACAGGAGAACGGCCTTACCGCTGTTCCCTCTGCCACGGGCTTTTAAGGCCTTG
GCACCTCTGCTCCGGCACCAACACCGACATGGGGTAGAGCCCGGCACCTCTGAAAGGCTTCTACCCACAA
CAACAACCGGACAGCCAACTCAAGGGTGGCCAGGAGAGGTGGAGGTGGTGTGGCTGCAGCAGCTGC
GGGTGCTGTGGTAGGCAAGCCTTTCGCCTGCAGGTTTTGTGCCAAGCCCTCCGACGCTCCTCAGATATG
CGAGACCAGAGCGGGTGCATACAGGGGAGCGGCCCTACCACTGTAGTATCTGTGGCAAAGGCTTACAGC
AGTCCTCTGTGCTGAGTGGTCATGCCCGCATCCACACAGGCGAGCGCCCTTCCGATGCATGCTCTGTGA
CCGCACCTTCAATAACTCCTCCAACCTCCGCAACACAGCGCACCCACTTCCATGGACCAGGGTCTGGA
GTGGGAGAGTCTAGAGGCCAACTGAGATCATCTTCAGTGTCCCAAGAATCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR225011 representing NM_001039532 Red=Cloning site Green=Tags(s) MAAAPDPPIPSSPTRESPTSPEPPDLVLPDGRPVTPPGDLIEIQVVKVTDTTSVPEPPEPGSFHCALCP AAFRLVSELLFHEHGLASMEGLGQGDPSRCHVCGHSCPGPASLRAHYSLHTGERPYRCSLCPRAFKAL APLLRHQHRHGVPGTSERLPTTTTGQPNRVAQERSEVVMAAAAGAVVGKPFACRFCAKPFRRSSDM RDHERVHTGERPYHCSICGKFTQSSVLSGHARIHTGERPFRCLCDRTFNSSNFRKHQRTHFHGPGSG VGESRGQLRSSVSQES TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_001039532
ORF Size:	891 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001039532.2, NP_001034621.1</u>
RefSeq Size:	2146 bp
RefSeq ORF:	894 bp
Locus ID:	654801
UniProt ID:	<u>Q8BI69</u>
Cytogenetics:	7 A1
MW:	32.2 kDa