

Product datasheet for **MR225016L4V**

Rfx3 (NM_001166414) Mouse Tagged ORF Clone Lentiviral Particle

Product data:


Product Type:	Lentiviral Particles
Product Name:	Rfx3 (NM_001166414) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rfx3
Synonyms:	C230093O12Rik; MRFX3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001166414
ORF Size:	2250 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225016).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001166414.1
RefSeq Size:	9146 bp
RefSeq ORF:	2250 bp
Locus ID:	19726
UniProt ID:	P48381
Cytogenetics:	19 C1



[View online »](#)

Gene Summary:

Transcription factor required for ciliogenesis and islet cell differentiation during endocrine pancreas development. Essential for the differentiation of nodal monocilia and left-right asymmetry specification during embryogenesis. Required for the biogenesis of motile cilia by governing growth and beating efficiency of motile cells (PubMed:15121860, PubMed:19671664). Also required for ciliated ependymal cell differentiation (PubMed:16930429). Together with RFX6, participates in the differentiation of 4 of the 5 islet cell types during endocrine pancreas development, with the exception of pancreatic PP (polypeptide-producing) cells (PubMed:17229940). Regulates transcription by forming a heterodimer with another RFX protein and binding to the X-box in the promoter of target genes (By similarity). Regulates the expression of genes involved in ciliary assembly (DYNC2LI1, FOXJ1 and BBS4) and genes involved in ciliary motility (DNAH11, DNAH9 and DNAH5). Represses transcription of MAP1A in non-neuronal cells but not in neuronal cells. [UniProtKB/Swiss-Prot Function]

Product images:

[MR225016L4] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with MR225016L4V particle to overexpress human Rfx3-mGFP fusion protein.