

## Product datasheet for **MR225021L3V**

### Wiz (NM\_011717) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Wiz (NM_011717) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Wiz
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_011717
ORF Size:	2865 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225021).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_011717.4</a>
RefSeq Size:	4247 bp
RefSeq ORF:	2868 bp
Locus ID:	22404
UniProt ID:	<a href="#">O88286</a>
Cytogenetics:	17 B1
Gene Summary:	May link EHMT1 and EHMT2 histone methyltransferases to the CTBP corepressor machinery. May be involved in EHMT1-EHMT2 heterodimer formation and stabilization. [UniProtKB/Swiss-Prot Function]



[View online »](#)