

## Product datasheet for **MR225025L4V**

### Nasp (NM\_016777) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Nasp (NM_016777) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nasp
Synonyms:	5033430J04Rik; AI131596; AI317140; D4Ertd767e; Epcs32; Nasp-T
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_016777
ORF Size:	2319 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225025).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_016777.3</a> , <a href="#">NP_058057.3</a>
RefSeq Size:	3087 bp
RefSeq ORF:	2322 bp
Locus ID:	50927
Cytogenetics:	4 53.24 cM
Gene Summary:	Required for DNA replication, normal cell cycle progression and cell proliferation. Forms a cytoplasmic complex with HSP90 and linker H1 histones and stimulates HSP90 ATPase activity. NASP and H1 histone are subsequently released from the complex and translocate to the nucleus where the histone is released for binding to DNA.[UniProtKB/Swiss-Prot Function]



[View online »](#)