

## Product datasheet for **MR225061L4V**

### Eri1 (NM\_026067) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Eri1 (NM_026067) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Eri1
Synonyms:	3'hexo; 3110010F15Rik; eri-1; Thex1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_026067
ORF Size:	1035 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225061).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_026067.3</a> , <a href="#">NP_080343.4</a>
RefSeq Size:	5067 bp
RefSeq ORF:	1038 bp
Locus ID:	67276
UniProt ID:	<a href="#">Q7TMF2</a>
Cytogenetics:	8 A4



[View online »](#)

**Gene Summary:**

RNA exonuclease that binds to the 3'-end of histone mRNAs and degrades them, suggesting that it plays an essential role in histone mRNA decay after replication. A 2' and 3'-hydroxyl groups at the last nucleotide of the histone 3'-end is required for efficient degradation of RNA substrates. Also able to degrade the 3'-overhangs of short interfering RNAs (siRNAs) in vitro, suggesting a possible role as regulator of RNA interference (RNAi). Binds with high affinity to the 3' side of the stem-loop structure and to the downstream cleavage product (DCP) of histone pre-mRNAs. Requires for binding the 5'-ACCCA-3' sequence present in stem-loop structure. Able to bind other mRNAs (By similarity). Required for 5.8S rRNA 3'-end processing. Also binds to 5.8s ribosomal RNA.[UniProtKB/Swiss-Prot Function]