

Product datasheet for **MR225065L4V**

Slc25a13 (NM_001177572) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc25a13 (NM_001177572) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Slc25a13
Synonyms:	AI785475; Ctrn
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001177572
ORF Size:	2025 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225065).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001177572.1 , NP_001171043.1
RefSeq Size:	3130 bp
RefSeq ORF:	2028 bp
Locus ID:	50799
UniProt ID:	Q9QXX4
Cytogenetics:	6 A1
Gene Summary:	Mitochondrial and calcium-binding carrier that catalyzes the calcium-dependent exchange of cytoplasmic glutamate with mitochondrial aspartate across the mitochondrial inner membrane. May have a function in the urea cycle.[UniProtKB/Swiss-Prot Function]



[View online »](#)