

Product datasheet for **MR225088**

Fnbp1 (NM_001177648) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fnbp1 (NM_001177648) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Fnbp1
Synonyms:	1110057E06Rik; 2210010H06Rik; FBP1; Fbp17
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>MR225088 representing NM_001177648
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAGCTGGGCACTGAGCTCTGGGATCAGTTTGACAACCTGGAAAAACATACACAGTGGGGAATCGATA
 TTCTTGAGAAATACATCAAGTTTGTCAAGGAGAGGACGGAGATTGAGCTCAGCTATGCCAAGCAACTCAG
 GAATCTTTCAAAGAAATACCAACCTAAGAAGAACTCGAAGGAAGAGGAGGAGTACAAGTACACGGCTTGC
 AAAGCCTTTCTTCCACCTGAATGAGATGAATGACTACGCCGGGCAGCAGAGGTCATCTCTGAGAACA
 TGACGTCACAGATCACGGTGGACCTGATGCGCTACGTTACAGGAGCTGAAGCAGGAGAGGAAATCGAACTT
 CCATGATGGACGGAAGGCTCAGCAGCACATAGAAACGTGTTGGAAGCAACTGGAGTCAAGTAAAGAGGAGG
 TTTGAGCGGGACTGTAAAGAACCCGACCGGCACAGCAGTACTTCGAGAAAATGGACGCTGACATCAACG
 TGACCAAGCGGATGTGAAAAGGCACGACAACAAGCTCAGATACGCCAGCAAATGGCAGAGGACAGCAA
 AGCAGATTACTCCTTAATCCTGCAGAGATCAACCAGGAGCAATGGGAATACTACCATACCCACATCCCC
 AACATCTTCCAGAAAATACAAGAGATGGAGGAGAGGCGGATTGTGCGTATTGGGGAGTCCATGAAGACGT
 ACGCAGAGGTGGACCGGACGGTATACCCATCATCGGGAAGTGCCTGGACGGGATAGTGAAGCGGCCGA
 GTCTATCGACCAGAAAAACGACTCCCAGCTGGTCTAGAACCTATAAGTCAGGATTCGAGCCTCCTGGA
 GACATTGAATTCGAAGATTACACACAGCCAATGAAACGCACAGTGTGAGCAACAGCCTTTCCAGCTCCA
 AAGAAGGCAAGCTGAGCTCAGATTTGGCGCAAGTCCAGAGGCAAGCTCTGGCCATTATCAAGAAAAA
 CAAGCTTATGTCCTTTTAAATCCCCCATCAGCCTCCCCCTCCGCCCTGCCTGCCTCACCTCT
 GCTGTTCCCAACGGCCCCAGTCTCCAAGCAGCAAAGGAACCCCTCTCCACCGCTTCAACGAGTTCA
 TGACCTCAAACCCAAAATCCACTGCTTCCGGAGCCTAAAGCGTGGGCTTCTCTCAAGCTGGGTGTAC
 ACCAGAAGACTTCAGCAACTTCCACCTGAGCAGAGAAAGGAAAAAATAACAGAAAAGTTGACGATCTC
 AATAGAGAGATACAGAAGGAGACGGACAGAGAGATGCCATACCAAAAATGAAAGATGTGTACCTAAAGA
 ACCCTCAGATGGGAGACCCAGCCAGCCTGGACCAGAAGCTCACCGAAGTCACCCAGAACATAGAGAACT
 GCGGCTGGAGGCTCAGAAGTTTGGGCTGGCTGGCTGAGGTAGAAGGCAGACTCCAGCTCGGAGTGAG
 CAGGCACGCCGAGAGTGGACTGTATGATGGCCAGACACACCAGACGGTCACTAACTGTGCCAGGACC
 GGGAGAGCCAGATGGTAGTTACACAGAGGAGCAAAGCCAGGAGAGCGAGCACAAGGTCCTGGCCCCGA
 TTTGACGATGAATTTGATGATGAGGAGCCGCTCCAGCCATAGGACCTGCAAGGCCCTACACATTT
 GAAGGTCAGAACGAAGGCCACATTTCCGTAGTTGAAGGAGAGACGCTGAGCGTATTGAAGAGGACAAAG
 GCGATGGGTGACTCGCATCCGCAGAAAATGAAGACGAGGAGGTTACGTCCCACCTCTACGTCGAAGT
 CTATTTAGACAAAACGCCAAAGTTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>MR225088 representing NM_001177648
 Red=Cloning site Green=Tags(s)

MSWGTTELWDQFDNLEKHTQWIDILEKYIKFVKERTEIELSYAKQLRNL SKKYQPKKNSKEEEYKYTAC
 KAFLSTLNEMNDYAGQHEVISENMTSQITVDMRYVQELKQERKSNFHDGRKAQQHIETCWKQLESSKRR
 FERDCKEADRAQQYFEKMDADINVTKADVEKARQQAQIRQQMAEDSKADYSLILQRFNQEQWEYHHTHIP
 NIFQKIQEMEERRIVRIGESMKTYAEVDRQVPIPIIGKCLDGIKAAESIDQKNSQLVVEAYKSGFEPFG
 DIEFEDYTQPMKRTVSDNSLSSSKEGKPELRFGGKSRGKLWPF IKKNKLSLLTSPHQPPPPPPASASPS
 AVPNGPQSPKQKPEPLSHRFNEFMTSKPKIHCFRSLKRGLSLKLGVTPEDFSNFPPEQRRKLLQKVVDDL
 NREIQKETDQRDAITKMKDVYLNPKQMGDPASLDQKLTEVTQNIKLRLEAQKFEAWLAEVEGRLPARSE
 QARRQSGLYDGQTHQVTNCAQDRESPDGSYTEEQSQSEHKVLPDFDDEFDDEEPLPAIGTKALYTF
 EQNQEGLISVVEGETLSVIEEDKGDGWTRIRRNEDEEGYVPTSIVEVYLDKNAKGS


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001177648

ORF Size: 1848 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001177648.1 , NP_001171119.1
RefSeq Size:	4743 bp
RefSeq ORF:	1851 bp
Locus ID:	14269
UniProt ID:	Q80TY0
Cytogenetics:	2 B
MW:	71.8 kDa
Gene Summary:	Required to coordinate membrane tubulation with reorganization of the actin cytoskeleton during the late stage of clathrin-mediated endocytosis. Binds to lipids such as phosphatidylinositol 4,5-bisphosphate and phosphatidylserine and promotes membrane invagination and the formation of tubules. Also enhances actin polymerization via the recruitment of WASL/N-WASP, which in turn activates the Arp2/3 complex. Actin polymerization may promote the fission of membrane tubules to form endocytic vesicles. May act as a link between RND2 signaling and regulation of the actin cytoskeleton. May be required for the lysosomal retention of FASLG/FASL (By similarity).[UniProtKB/Swiss-Prot Function]