

Product datasheet for **MR225121**

Fgfr2 (NM_201601) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fgfr2 (NM_201601) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Fgfr2
Synonyms:	AU043015; AW556123; Bek; Fgfr-2; Fgfr-7; Fgfr7; KGFR; KGFRTr; sv5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR225121 representing NM_201601
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGGATTACCGTCCACGTGGAGATATGGAAGAGGACCAGGGATTGGCACTGTGACCATGGTCAGCTGGG
 GCGCTTCATCTGCCTGGTCTTGGTCACCATGGCAACCTTGTCCCTGGCCCGGCCCTCCTCAGTTTAGT
 TGAGGATACCACTTTAGAACAGAGGAGCACCCTACTGGACCAACACCGAGAAGATGGAGAAGCGGCTC
 CACGCTGTCCCTGCCCAACACTGTGAAGTCCGCTGTCCGGCTGGGGGAATCCAACGCCACAATGA
 GGTGGTTAAAAACGGGAAGGAGTTTAAGCAGGAGCATCGCATTGGAGGCTATAAGGTACGAAACCAGCA
 CTGGAGCCTTATTATGAAAAGTGTGGTCCCGTCAGACAAAGGCAACTACACCTGCCTGGTGGAGAATGAA
 TACGGGTCCATCAACCACACCTACCACCTCGATGTCGTTGAACGGTCACCACACCGGCCCATCTCCAAG
 CTGGACTGCCTGCAAATGCCTCCACGGTGGTGGAGGGGATGTGGAGTTTGTCTGCAAGGTTTACAGCGA
 TGCCACAGCCACATCCAGTGGATCAAGCACGTGGAAAAGAACGGCAGTAAATACGGGCTGATGGGCTG
 CCCTACCTCAAGTCTGAAGCACTCGGGGATAAATAGCTCCAATGCAGAAGTGCTGGCTCTGTTCAATG
 TGACGGAGATGGATGCTGGGGAATATATGTAAGGTCTCCAATTATATAGGGCAGGCCAACAGTCTGC
 CTGGCTCACTGTCTGCCAAACAGCAAGCGCCTGTGAGAGAGAAGGAGATCACGGCTTCCCCAGATTAT
 CTGGAGATAGCTATTTACTGCATAGGGGTCTTCTTAATCGCCTGCATGGTGGTGACAGTCATCTTTTGCC
 GAATGAAGACCACGACCAAGAAGCCAGACTTCAGCAGCCAGCCAGCTGTGCACAAGCTGACCAAGCGCAT
 CCCCTGCGGAGACAGGTAACAGTTTCGGCCGAGTCCAGCTCCTCCATGAACTCCAACACCCCGTGGTG
 AGGATAACAACCGCTCTGTCCTCAACAGCGGACACCCCGATGCTAGCAGGGGTCTCCGAGTATGAGTTGC
 CAGAGGATCCAAAGTGGGAATTCACAGAGATAAGCTGACGCTGGGCAAACCCCTGGGGGAAGTTGCTT
 CGGGCAAGTAGTCATGGCTGAAGCAGTGGGAATCGATAAAGACAAACCAAGGAGCGGTACCCGTGGCA
 GTGAAGATGTTGAAAGATGATGCCACAGAGAAGGACCTGTCTGATCTGGTATCAGAGATGGAGATGATGA
 AGATGATTGGAAACATAAGAACATTATCAACCTCCTGGGGCCTGCACGCAGGATGGACCTCTCTACGT
 CATAGTTGAATATGCATCGAAAGGCAACCTCCGGGAATACCTCCGAGCCCGGAGGCCACCTGGCATGGAG
 TACTCCTATGACATTAACCGTGTCCCGAGGAGCAGATGACCTTCAAGGACTTGGTGTCTGCACCTACC
 AGCTGGCTAGAGGCATGGAGTACTTGGCTTCCAAAAATGTATCCATCGAGATTTGGCTGCCAGAAACGT
 GTTGGTAACAGAAAACAATGTGATGAAGATAGCAGACTTTGGCCTGGCCAGGGATATCAACAACATAGAC
 TACTATAAAAAGACCACAAATGGGCGACTTCCAGTCAAGTGGATGGCTCCTGAAGCCCTTTTGGATAGAG
 TTTACACTCATCAGAGCGATGTCTGGTCCCTCCGGGTGTTAATGTGGGAGATCTTACTTTAGGGGGCTC
 ACCCTACCCAGGGATTCCCGTGGAGGAACTTTTAAGCTGCTCAAAGAGGGACACAGGATGGACAAGCCC
 ACCAACTGCACCAATGAACTGTACATGATGATGAGGGATTGCTGGCATGCTGTACCTCAGAGACCCCA
 CATTCAAGCAGTTGGTCGAAGACTTGGATCGAATTCTGACTCTCACAAACCAATGAGGAATACTTGGATCT
 CACCCAGCCTCTCGAACAGTATTCTCCTAGTTACCCCGACACAAGGAGCTCTTGTCTTCAGGGGACGAT
 TCTGTGTTTTCTCCAGACCCCATGCCTTATGAACCCTGTCTGCCTCAGTATCCACACATAAACGGCAGTG
 TAAAAACA

AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR225121 representing NM_201601
 Red=Cloning site Green=Tags(s)

MGLPSTWRYGRGPGIGTVMVSWGRFICLVLTMTLSLARPFSLVEDTTLEPEGAPYWTNTEKMEKRL
 HAVPAANTVKFRCPAGGNPTPTMRWLKNGKEFKQEHRIGGYKVRNQHWSLIMESVVP SDKGN YTC L VENE
 YGSINHTYHLDVVERSPHRPILQAGLPANASTVVGDDVEFVCKVYSDAQPHIQWIKHVEKNGSKYGPDGL
 PYLKVLKHSGINSNAEVLALFNVTENDAGEYICKVSNYIGQANQSAWLTVLPKQQAPVREKEITASPDY
 LEIAIYCIQVFLIACMVVTVIFCRMKTTTKKPDFSSQPAVHKLTKRIPLRRQVTVSAESSMSNTPLV
 RITTRLSSTADTPMLAGVSEYELPEDPKWEFPRDKLTLGKPLGEGCFQVVMMAEAVGIDKDKPEAVTVA
 VKMLKDDATEKDLSDLVSEMEMMKMIGKHKNIINLLGACTQDGPLYVIVEYASKGNLREYLRRRPPGME
 YSYDINRVPEEQMTFKDLVSCTYQLARGMEYLASQKCIHRDLAARNVLTENNVMKIADFLARDINNID
 YYKKTNGRLPVKWMPEALFDRVYTHQSDVWSFGVLMWEIFTLGGSPYPGIPVEELFKLLKEGHRMDKP
 TNCTNELYMMMRDCWHAVPSQRP TFKQL VEDLDRILTLTNEEYLDLTQPLEQYSPSPDTRSSCSGDD
 SVFSPDPMPYEPCLPQYPHINGSVKT



SGP TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9030_e10.zip

Restriction Sites: SgfI-RsrII

Cloning Scheme:

Cloning sites used for ORF Shutting:



* The last codon before the Stop codon of the ORF

ACCN: NM_201601


ORF Size: 2178 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_201601.2</u> , <u>NP_963895.2</u>
RefSeq Size:	4881 bp
RefSeq ORF:	2181 bp
Locus ID:	14183
UniProt ID:	<u>P21803</u>
Cytogenetics:	7 F3
MW:	82.1 kDa
Gene Summary:	<p>Tyrosine-protein kinase that acts as cell-surface receptor for fibroblast growth factors and plays an essential role in the regulation of cell proliferation, differentiation, migration and apoptosis, and in the regulation of embryonic development. Required for normal embryonic patterning, trophoblast function, limb bud development, lung morphogenesis, osteogenesis and skin development. Plays an essential role in the regulation of osteoblast differentiation, proliferation and apoptosis, and is required for normal skeleton development. Promotes cell proliferation in keratinocytes and immature osteoblasts, but promotes apoptosis in differentiated osteoblasts. Phosphorylates PLCG1, FRS2 and PAK4. Ligand binding leads to the activation of several signaling cascades. Activation of PLCG1 leads to the production of the cellular signaling molecules diacylglycerol and inositol 1,4,5-trisphosphate. Phosphorylation of FRS2 triggers recruitment of GRB2, GAB1, PIK3R1 and SOS1, and mediates activation of RAS, MAPK1/ERK2, MAPK3/ERK1 and the MAP kinase signaling pathway, as well as of the AKT1 signaling pathway. FGFR2 signaling is down-regulated by ubiquitination, internalization and degradation. Mutations that lead to constitutive kinase activation or impair normal FGFR2 maturation, internalization and degradation lead to aberrant signaling. Over-expressed FGFR2 promotes activation of STAT1.[UniProtKB/Swiss-Prot Function]</p>

Product images:



Circular map for MR225121