

## Product datasheet for **MR225130L3V**

### Krtdap (NM\_001033131) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Krtdap (NM_001033131) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Krtdap
Synonyms:	Kdap; sk221
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001033131
ORF Size:	306 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225130).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001033131.3</a>
RefSeq Size:	456 bp
RefSeq ORF:	309 bp
Locus ID:	64661
UniProt ID:	<a href="#">Q3V2T4</a>
Cytogenetics:	7 B1
Gene Summary:	May act as a soluble regulator of keratinocyte differentiation. May play an important role in embryonic skin morphogenesis (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)