

## Product datasheet for **MR225142L4V**

### **Dpp6 (NM\_010075) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Dpp6 (NM_010075) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dpp6
Synonyms:	B930011P16Rik; D5Buc3; D5Buc4; D5Buc5; Dpp-6; DPP VI; Gm1377; In(5)6H-p; Peplb; Rw
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_010075
ORF Size:	2412 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225142).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_010075.2</a> , <a href="#">NP_034205.1</a>
RefSeq Size:	4866 bp
RefSeq ORF:	2415 bp
Locus ID:	13483
UniProt ID:	<a href="#">Q9Z218</a>
Cytogenetics:	5 12.92 cM



[View online »](#)

**Gene Summary:**

Promotes cell surface expression of the potassium channel KCND2 (PubMed:22311982). Modulates the activity and gating characteristics of the potassium channel KCND2 (PubMed:22311982). Has no dipeptidyl aminopeptidase activity (By similarity). [UniProtKB/Swiss-Prot Function]