

Product datasheet for MR225149L4

Sox15 (NM_009235) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Sox15 (NM_009235) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Sox15
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225149).
Restriction Sites:	Sgfl-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_009235
ORF Size:	693 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_009235.2 , NP_033261.1
RefSeq Size:	1156 bp
RefSeq ORF:	696 bp
Locus ID:	20670
UniProt ID:	P43267
Cytogenetics:	11 42.86 cM
Gene Summary:	Transcription factor that binds to DNA at the 5'-ACAATG-3' consensus sequence (PubMed:10821863, PubMed:15863505, PubMed:16759287, PubMed:17363903). Acts as a transcriptional activator and repressor (PubMed:10821863, PubMed:15863505, PubMed:16759287). Binds synergistically with POU5F1 (OCT3/4) to gene promoters (PubMed:15863505). Binds to the FOXC1 promoter and recruits FHL3, resulting in transcriptional activation of FOXC1 which leads to myoblast proliferation (PubMed:17363903). Acts as an inhibitor of myoblast differentiation via transcriptional repression which leads to down-regulation of the muscle-specific genes MYOD and MYOG (PubMed:10821863). Involved in trophoblast giant cell differentiation via enhancement of HAND1 transcriptional activity (PubMed:16759287). Regulates transcription of HRC via binding to its proximal enhancer region (PubMed:15863505). Involved in skeletal muscle regeneration (PubMed:15367664, PubMed:17363903). Also plays a role in the development of myogenic precursor cells (PubMed:15367664).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR225149L4