

Product datasheet for **MR225168L3V**

Mkx (NM_177595) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mkx (NM_177595) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mkx
Synonyms:	9430023B20Rik; Irx11
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_177595
ORF Size:	1065 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225168).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_177595.4 , NP_808263.2
RefSeq Size:	3201 bp
RefSeq ORF:	1065 bp
Locus ID:	210719
UniProt ID:	Q8BIA3
Cytogenetics:	18 4.53 cM
Gene Summary:	May act as a morphogenetic regulator of cell adhesion. Participates in the early events that lead to differentiation.[UniProtKB/Swiss-Prot Function]



[View online »](#)