

## Product datasheet for **MR225234L3V**

### Cldn24 (NM\_001111318) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Cldn24 (NM_001111318) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cldn24
Synonyms:	Gm10107
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001111318
ORF Size:	660 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225234).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001111318.1</a> , <a href="#">NP_001104788.1</a>
RefSeq Size:	663 bp
RefSeq ORF:	663 bp
Locus ID:	100039801
Cytogenetics:	8 B1.1



[View online »](#)

**Gene Summary:**

This gene encodes a member of the claudin family. Claudins are integral membrane proteins and components of tight junction strands. Tight junction strands serve as a physical barrier to prevent solutes and water from passing freely through the paracellular space between epithelial or endothelial cell sheets, and also play critical roles in maintaining cell polarity and signal transductions. The protein encoded by this gene is 75% identical to the human homolog. Similar to the human gene, this gene is upstream of the Cldn22 gene, which overlaps the Wwc2 gene on the opposite strand. [provided by RefSeq, Aug 2010]