

## Product datasheet for **MR225266L3V**

### **Pou3f2 (NM\_008899) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Pou3f2 (NM_008899) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pou3f2
Synonyms:	9430075J19Rik; A230098E07Rik; Brn-2; Brn2; oct-7; OTF-7; Otf7
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008899
ORF Size:	1335 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225266).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_008899.1</a> , <a href="#">NP_032925.1</a>
RefSeq Size:	1338 bp
RefSeq ORF:	1338 bp
Locus ID:	18992
UniProt ID:	<a href="#">P31360</a>
Cytogenetics:	4 9.73 cM



[View online »](#)

**Gene Summary:**

Transcription factor that plays a key role in neuronal differentiation (PubMed:24243019). Binds preferentially to the recognition sequence which consists of two distinct half-sites, ('GCAT') and ('TAAT'), separated by a non-conserved spacer region of 0, 2, or 3 nucleotides (By similarity). The combination of three transcription factors, ASCL1, POU3F2/BRN2 and MYT1L, is sufficient to reprogram fibroblasts and other somatic cells into induced neuronal (iN) cells in vitro (PubMed:20107439, PubMed:24243019, PubMed:27281220). Acts downstream of ASCL1, accessing chromatin that has been opened by ASCL1, and promotes transcription of neuronal genes (PubMed:24243019).[UniProtKB/Swiss-Prot Function]