

Product datasheet for **MR225267L3V**

Mcoln2 (NM_026656) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mcoln2 (NM_026656) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mcoln2
Synonyms:	3300002C04Rik; AI549968; C86638; TRPML2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_026656
ORF Size:	1698 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225267).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026656.4
RefSeq Size:	2519 bp
RefSeq ORF:	1701 bp
Locus ID:	68279
UniProt ID:	Q8K595
Cytogenetics:	3 H2



[View online »](#)

Gene Summary:

Nonselective cation channel probably playing a role in the regulation of membrane trafficking events. Acts as Ca(2+)-permeable cation channel with inwardly rectifying activity (PubMed:19763610). May activate ARF6 and be involved in the trafficking of GPI-anchored cargo proteins to the cell surface via the ARF6-regulated recycling pathway (By similarity). May play a role in immune processes. In adaptive immunity, TRPML2 and TRPML1 may play redundant roles in the function of the specialized lysosomes of B cells (PubMed:17050035). In the innate immune response, may play a role in the regulation of chemokine secretion and macrophage migration (PubMed:26432893). Through a possible and probably tissue-specific heteromerization with MCOLN1 may be at least in part involved in many lysosome-dependent cellular events.[UniProtKB/Swiss-Prot Function]