

Product datasheet for **MR225297L4V**

Ly96 (NM_016923) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ly96 (NM_016923) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ly96
Synonyms:	ESOP-1; MD-2; MD2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_016923
ORF Size:	483 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225297).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016923.2 , NP_058619.1
RefSeq Size:	603 bp
RefSeq ORF:	483 bp
Locus ID:	17087
UniProt ID:	Q9JHF9
Cytogenetics:	1 A3



[View online »](#)

Gene Summary:

Binds bacterial lipopolysaccharide (LPS) (PubMed:22532668). Cooperates with TLR4 in the innate immune response to bacterial lipopolysaccharide (LPS), and with TLR2 in the response to cell wall components from Gram-positive and Gram-negative bacteria. Enhances TLR4-dependent activation of NF-kappa-B. Cells expressing both LY96 and TLR4, but not TLR4 alone, respond to LPS (PubMed:10725698).[UniProtKB/Swiss-Prot Function]