

Product datasheet for **MR225356L4V**

Pla2g1b (NM_011107) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pla2g1b (NM_011107) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pla2g1b
Synonyms:	Pla2a; sPLA2IB
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_011107
ORF Size:	441 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR225356).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011107.1
RefSeq Size:	552 bp
RefSeq ORF:	441 bp
Locus ID:	18778
UniProt ID:	Q9Z0Y2
Cytogenetics:	5 F
Gene Summary:	PA2 catalyzes the calcium-dependent hydrolysis of the 2-acyl groups in 3-sn-phosphoglycerides, this releases glycerophospholipids and arachidonic acid that serve as the precursors of signal molecules.[UniProtKB/Swiss-Prot Function]


[View online »](#)